

U0099.S – ELECTRIC BED



GENERAL SPECIFICATION

- Electric bed with four functions.
- Dimensions: 2120 / 970 / 500 – 820 mm.
- Maximum load: 120kg.
- Four-section.
- Steel frame and bases of individual sections made of galvanized steel mesh - powder coating.
- The head and foot boards are plastic with decoration plate.
- Movement of sections via remote control:
 - back lift 0°-75°.
 - lower leg lift 0°-40°.
 - height lift - 320 mm.
 - base tilt: trendelenburg/anti-trendelenburg.
- Electric motors - 4 pcs. Electronic unit - 1 pc.
- Remote control:
 - 2 pcs. buttons for adjusting the back section.
 - 2 pcs. buttons for adjusting sections at the lower leg.
 - 2pcs. buttons for simultaneous movement of the sections.
 - 2pcs. buttons for height adjustment.
 - 2pcs. buttons for trendelenburg/anti-trendelenburg.
- *Aluminum side panels with a mechanism for removing and lifting.
- 4pcs. wheels rotating 360° all with sperate foot brake.

IHC

MATTRESS SPECIFICATION

Dimensions:2000x880x120mm.

- foam

- anti-decubitus. 5 ZONE.

Cover fabric:

- Double elastic polyurethane coating;

- The fabric meets the highest requirements for optimal weight and resistance to the smallest particles;

- Waterproof and "breathable" fabric (vapor permeable);

- Antibacterial, hypoallergenic;

- Washing at a temperature of 95°C;

- Withstands modern disinfectants.

- With zipper (on the short and long side).

- Eco-friendly.

- CRIB 5 fireproof.



IHC Medikal A.S.

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

www.ihcmed.com.tr

U0013.S – MECHANICAL HOSPITAL BED



TECHNICAL SPECIFICATION

- Mechanical bed with two functions.
- Dimensions: 2120/970/500mm.
- Maximum safe load: 120 kg.
- Four-section inner frame.
- Steel supporting frame and bases of individual sections made of galvanized steel mesh and metal profiles - powder coating.
- ABS trays with two handles - removable, with the possibility of placing a plate. Plastic head and foot boards.
- Lifting of the sections using 2 pcs. manual screw mechanisms.
- Movement of the sections:
 - lifting at the back 0°-75°.
 - lifting at the lower leg/calf 0°-40°.
- 4 pcs. wheels ϕ 125mm, rotating 360° each with foot brake – isolated with a plastic fork preventing contamination.
- Side barriers with a mechanism for removing and lifting.
- Adjustable / in height / stainless steel infusion stand with the possibility of placing on both sides of the bed.
- Without mattress.



U0026.S – HOSPITAL CABINETS WITH TABLE



TECHNICAL SPECIFICATION

- The dimensions are 510 mm × 390 mm × 770 mm.
- The material is metal and chipboard.
- The table is with mechanism.
- It is on wheels.
- There is a drawer.
- There is a cabinet.

