

To: MBAL Sveti Ivan Rilski


Dupnitsa, Republic of Bulgaria



Re: announcement No. 00565-2023-0002

for public procurement procedure

DECLARATION

We **IHC Medikal A.S.**, with address Bahcelievler District Yavuz Road No:31 Hasanoglan Elmadag, Ankara Turkey , hereby declare and confirm that our products, offered through our distributor **ELPAK Leasing EOOD** with business address 12 Ivan Bogorov str., 9002 Varna, Bulgaria, in Public Tender (public procurement procedure) as per announcement No. **00565-2023-0002** by **MBAL Sveti Ivan Rilski, Dupnitsa, Bulgaria (as 'contracting authority')** fully covers the procedure's requirements, the technical parameters and will be designed especially for the tender with the specification as follows:

PRODUCT CODE	PRODUCT PHOTO	PRODUCT NAME	DESCRIPTION	QUAN
U0014.T4.B2. K0.YWM.WO U0018.S		HOSPITAL BED WITH 1 CRANK MATTRESS	No side boards. Back section adjustment: 0°-70° (by crank). Powder coated construction. Anti-shock elements in the four corners of the bed. Wheels, brakes. Maximum load capacity: 170 kg. Overall dimensions: 2100 x 980 x h 500 mm. Mattress: Waterproof, antibacterial case, blue colour, with a zipper. Filling: foam. Dimensions: 1900 x 900 x 140 mm.	46 Pc

<p>U0003.S.B5.T 16S.K5.E0.M1 .M3.C1.YWM. W1</p> <p>U0018</p>		<p>HOSPITAL BED 4 MOTORS</p> <p>MATTRESS</p>	<p>Electric resuscitation bed with mattress: Steel frame and bases of the individual sections of steel galvanized mesh - powder coating. ABS boards with two handles - removable. Movement of sections by remote control: Elevation at the back 0°-85° Lower leg elevation 0°-45°. Calf raise by mechanical. Lifting height – 320 mm. Base slope: trendelenburg/ reverse trendelenburg. Electric motors - 4 pieces; Electronic block - 1 piece. Remote control: - 2 buttons for adjusting the back section. - 2 buttons for adjusting sections at the lower leg. - 2 buttons for simultaneous movement of the sections. - 2 buttons for adjusting the height. - 2 trendelenburg/anti-trendelenburg buttons. ABS side partitions with lifting and lowering mechanism - 4 pieces. 4 pieces of 360° rotating wheels with foot brakes/4 pieces/. Adjustable/height/bedside stainless steel infusion stand. Bed dimensions: 2120 x 1000 x 520-840 mm.</p>	<p>4 Pc</p>
<p>U0024.S.T8.K L2</p>		<p>ABS BEDSIDE CABINET</p>	<p>4 wheels rotating 360°, with brakes on two of them. Two compartments. Drawer and cabinet. A drawer. Movable board from one side. Towel hanger - 4 pieces. Overall dimensions: 480x480x760 mm. 4 horizontal corners are aluminium profile. The other all parts are ABS plastic.</p>	<p>48 Pc</p>



We further declare and confirm that the devices shall fully comply with all of the above-mentioned technical specifications and features and all requirements of the contracting authority. The devices shall be designed and manufactured specifically for this public tender (public procurement procedure), for which the manufacturer possesses full certification documents, capacity and experience.

Please note that the product catalogues and other advertising materials of IHC Medikal A.S. contain only general presentation of our standard products, therefore the same should not be relied on by the contracting authority as regards the technical specification and features of the device subject to offer by our distributor ELPK Leasing EOOD.

IHC Medikal A.S.

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

www.ihcmed.com.tr