


To: MBAL Knyaginya Klementina

Sofia, Republic of Bulgaria

Re: announcement No. 00293-2023-0005

for public procurement procedure


CATALOG

PRODUCT CODE	PRODUCT PHOTO	PRODUCT NAME	DESCRIPTION	QUAN.
U0001.T24. K5.E1.M2. M3.C1.YMS. W1 U0018.S		HOSPITAL BED 4 MOTORS MATTRESS	Resuscitation hospital beds with mattress: 2 telescopic columns and possibility of fully electric bed drive. 4 sections - back, seat, thighs, legs. With electrically adjustable back, upper leg section, height, trendelenburg and reverse trendelenburg function. Back section lift from 0° to 70° (±2°). Leg section lift from 0° to 30° (±2°). Trendelenburg/Reverse Trendelenburg ±15°(±2°). Ability to raise the back and leg sections simultaneously. Adjustable bed height from 450 to 850 mm (±5 cm). Possibility of mechanical CPR for back rest movement. With 2 batteries. Metal mattress platform four sections, made of galvanised metal sheets, powder coated, with mattress ventilation openings and head and foot restraints; With removable plastic boards at the head and feet. Side rails - to consist of 2 x 2 medical HDPE sides, foldable and falling below the level of the mattress structure (in accordance	4 Pcs

IHC Medikal A.S.

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

www.ihcmed.com.tr

			<p>with European Standard EN60601-2-52), which can be lowered below the level of the mattress in its lowest position.</p> <p>The bed shall be protected by $\varnothing 75$ mm (± 3 mm) impact buffers located at the four corners of the bed.</p> <p>4 ABS wheels $\varnothing 150$ mm central braking castors one of them with directional lock.</p> <p>Dimensions of the mattress platform: 1970 x 860 mm (± 2 mm).</p> <p>External dimensions: 2180 x 1020 mm (± 2 mm);</p> <p>Loading capacity: 240 kg.</p> <p>Options to choose the head and foot boards and side rails laminated parts.</p> <p>Supplied with: anti-decubitus mattress with cover.</p> <p>The mattress of polyurethane foam, segmented laterally.</p> <p>Cover with zipper, made of waterproof fabric that is easy to clean and disinfection.</p>	
<p>U0099.S.T1 6.K1.YMS.W 0 U0018.S</p>		<p>HOSPITAL BED 4 MOTORS</p> <p>ANTI DECUBITUS MATTRESS</p>	<p>Resuscitation hospital beds with mattress: 4-section electric hospital bed with adjustable back rest, upper leg rest, height, trendelenburg and reverse trendelenburg by motors.</p> <p>With 4 electric motors to control the bed functions by remote control: Back section lift from 0° to 70° ($\pm 2^\circ$). Leg section lift from 0° to 35° ($\pm 2^\circ$).</p> <p>Option to raise the back section and the leg section simultaneously.</p> <p>Adjustable bed height from 420 to 820 mm (± 2 mm).</p> <p>Metal mattress platform with four sections, made of galvanised metal sheets, powder-coated, with mattress ventilation openings and head and foot restraints; The movable parts of the mattress platform shall allow easy cleaning;</p> <p>With removable impact resistant HDPE plastic boards at the head and foot.</p> <p>Side barriers - Metal tube side rails.</p> <p>With $\varnothing 75$ mm (± 3 mm) impact buffers located at all four corners of the bed.</p> <p>4 ABS wheels $\varnothing 125$ mm (± 2 mm), 2 with</p>	<p>2 Pcs</p>

IHC

			<p>brakes.</p> <p>Dimensions of the mattress platform: 2000 x 860 mm (± 2 mm).</p> <p>External dimensions: 2190 x 1005 mm (± 2 mm).</p> <p>Loading capacity: 230 kg.</p> <p>Option to choose of colour of head and foot boards laminates and the metal tube side rail RAL color.</p> <p>Supplied with: anti-decubitus mattress with cover.</p> <p>Mattress to be made of polyurethane foam, segmented laterally.</p> <p>Cover with zipper, made of waterproof fabric that is easy to clean and disinfection.</p>	
--	--	--	--	--