



| PRODUCT CODE | PRODUCT PHOTO | PRODUCT NAME | DESCRIPTION | QUAN. |
|---|---|--|---|---------------|
| <p>U0013.S. T36.B2.K 3S.YMS. W8S</p> <p>U0018.S</p> <p>U0112.T8 .KL2B (BLOW MOLDED)</p> |   | <p>MECHANICAL HOSPITAL BED</p> <p>MATTRESS</p> <p>HOSPITAL LOCKER</p> | <p>"Mechanical hospital bed, complete with a hospital mattress with a waterproof cover and a hospital locker."</p> <p>1. Mechanical bed with four sections, with a metal sheet and a surface with a special anti-corrosion coating.</p> <p>2. Head and foot boards are plastic with grip handles.</p> <p>3. Options to easily straighten the torso or legs in the lower leg and calf area of the patient.</p> <p>4. Be equipped with retractable sides to prevent the patient from falling. The sides equipped with locking mechanisms.</p> <p>5. 125 mm wheels equipped 2 with strong brakes.</p> <p>6. Features:</p> <ul style="list-style-type: none"> - Raising the patient's torso: from 0° to 75°. - Lower leg lift: from 0° to 40°. - Calf Raise. <p>7. To have an infusion stand, with 4 hooks.</p> <p>8. Patient weight 200 kg.</p> <p>9. Dimensions:</p> <ul style="list-style-type: none"> - Length: 2120 mm; ±20 mm. - Width: 900 mm; ±70 mm. - Height: 500 mm. <p>10. To be equipped with a medical mattress made entirely of foam with an easy-to-clean waterproof cover with a zipper.</p> <p>The mattress cover is fire resistant.</p> <p>11. Instructions for use in Bulgarian.</p> <p>12. Warranty period: 15 months including spare parts only.</p> <p>II. Hospital cabinet.</p> <ol style="list-style-type: none"> 1. Made of ABS plastic or equivalent, for easy cleaning and disinfection. 2. With a large cupboard below, with a shelf, a drawer above and a removable dining table/board. 3. There should be a hangers on the side. 4. It has 4 wheels - two of which have brakes. 5. Color - Grey. 6. Approximate Dimensions: 465 x 470 x 790 mm. ±5 mm. | <p>20 Pc.</p> |