




To: SBALAG Maichin Dom EAD

Sofia, Republic of Bulgaria




Re: announcement No. 00521-2023-0014

for public procurement procedure

## CATALOG

PRODUCT CODE	PRODUCT PHOTO	PRODUCT NAME	DESCRIPTION	QUAN.
NO 1 U0014.S. T4.B2.K0. YWM.W0		HOSPITAL BED WITH 1 CRANK	Two-section hospital bed: * To be with two-section, one mobile and one stationary; * To have steel frame and fundamentals of individual section of steel profile- powder painted construction; * To be able to raise the back section to 70°; * To be installation with soundproof wheels with rubberized coating and stop-system. * Overall dimensions up to 210/100/50 cm.	24 Pcs
NO 2 U0018.S		MATTRES S	* Anti-decubitus mattress, retrieved from high-quality- antibacterial, hypoallergenic, intended for hospital bed; * To be waterproof; * To be wearing case, with elastic polyurethane coating, for maintenance and cleaning with disinfectants; * Overall dimensions up to 200/90/15 cm.	24 Pcs
NO 3 U0024.S. T8.KL1		BEDSIDE CABINET	Attachment (hospital) locker: The device is made of material designed for hospital use and can with stand frequent treatment with disinfectants; * The device includes drawer, door * With wheels, with foot brakes- at two of them;	24 Pcs





# IHC

			* Size to 50/50/75 cm.	
<p><b>NO 4</b></p> <p><b>U0014.S.</b> <b>T4.B2.K0.</b> <b>YWM.WO</b></p>		<p><b>HOSPITAL BED WITH 1 CRANK</b></p>	<p>Two-section hospital bed:            * To be with two sections one mobile and one stationary;            * To have steel frame and fundamentals of individual sections of steel profile- powder painted construction;            * To be able to raise the back section to 70°;            * To be installation with soundproof wheels with rubberized coating and stop-system.            * Gauge size to 210/100/50 cm.</p>	<p><b>6 Pcs</b></p>
<p><b>NO 5</b></p> <p><b>U0099.S.</b> <b>T4S.B2.K0</b> <b>.M1.M3.C</b> <b>1.YWM.</b> <b>W8</b></p>		<p><b>HOSPITAL BED</b></p>	<p>Electric hospital bed with adjustable height:            Size: min 2100/950/500-820 mm.            * Steel frame and bases of the individual sections of galvanized steel mesh- powder painting.            * ABS trays with two handles, with the possibility of placing a plate.            * Lifting the sections by means of electric motor.            * Movement of sections by remote control:            - Lifting at the back 0°-75°.            - Lifting at the shin 0°-45°.            - Lifting at the height- 300 mm.            * Inclination of the base: trendelenburg/antitrendelenburg            Remote control with buttons for adjusting sections at the back, buttons for adjusting sections at the hip, buttons for simultaneous movement of sections, buttons for adjusting height and buttons for trendelenburg/antitrendelenburg.            * 4 wheels rotating on 360° with foot brakes- soundproofing.            * Stainless steel adjustable/height/infusion stand.            * Without side barriers.</p>	<p><b>2 Pcs</b></p>
<p><b>NO 6</b></p> <p><b>U0018.S</b></p>		<p><b>MATTRES S</b></p>	<p>* Antidekubital mattress made of high quality materials - antibacterial, hypoallergenic, designed for hospital bed;            * Be watertight;            * To be covered with a case, with a flexible polyurethane coating, for maintenance and cleaning with</p>	<p><b>8 Pcs</b></p>






**IHC Medikal A.S.**

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40


[www.ihcmed.com.tr](http://www.ihcmed.com.tr)

			<p>disinfectants;          * Overall dimensions up to 200/90/15 cm.</p>	
<p><b>NO 7</b>  <b>U0024.S</b>  <b>T8.KL1</b></p>		<p><b>BEDSIDE CABINET</b></p>	<p>Attachment (hospital) locker:          The device is made of material designed for hospital use and can with stand frequent treatment with disinfectants;          * The device includes drawer, door          * With wheels, with foot brakes- at two of them;          * Size to 50/50/75 cm.</p>	<p><b>8 Pcs</b></p>
<p><b>NO 8</b>  <b>U0047.S</b></p>		<p><b>MEDICAL CABINET</b></p>	<p>Double medical cabinet:          * To be used for storage of medical supplies;          * Cabinet to be of glazed doors for maximum visibility;          * Mechanical parts to ensure very cyclic operation;          * To be with adjustable legs, antibacterial base, neutral to disinfectants with easy to clean surface;          * The handles shall have a wear-resistant coating;          * With secure door or locking;          * Materials to withstand cleaning with disinfectants.          * Dimensions up to 90/45/ 180 cm.</p>	<p><b>2 Pcs</b></p>
<p><b>NO 9</b>  <b>U0049.S</b></p>		<p><b>CABINET</b></p>	<p>Double cabinet for medicines:          * The cabinet shall have an upper glazed part divided into three levels and a lower solid part with two levels;          * Mechanical parts to ensure very cyclic operation;          * To be with adjustable legs, antibacterial base, neutral to disinfectants with easy to clean surface;          * The handles shall have a wear-resistant coating;          * With secure door locking;          * Materials to withstand cleaning with disinfectants.          * Dimensions up to 90/45/ 180 cm.</p>	<p><b>1 Pcs</b></p>
<p><b>NO 10</b>  <b>U0114.S</b></p>		<p><b>GYNECOLOGICAL CHAIR</b></p>	<p>Gynecological electric chair:          * Electrostatic powder coating of components:          * Electric adjustment of the back section, height trendelenburg;          * Mattress platforms, made from</p>	<p><b>1 Pcs</b></p>

# IHC

			<p>polyurethane foam with high density;</p> <ul style="list-style-type: none"> <li>* Coating of artificial leather, sustainable of fire and disinfectants;</li> <li>* Kit: footrest, headrest, armrests, handles, waste container.</li> <li>* Workload: 180 kg.</li> </ul>	
NO 12 U0065.S		<b>EXAMINATION COUCH</b>	<p>Two-section medical couch:</p> <ul style="list-style-type: none"> <li>* Size 180/60/70cm.</li> <li>* Tubular construction of material allowing active disinfection;</li> <li>* To be upholstered in elastic foam and dressed in eco-leather to allow active disinfection</li> <li>* Be able to lift the back section.</li> <li>* Safe Workload: 120 kg.</li> </ul>	<b>2 Pcs</b>
NO 13 U0068.S		<b>SCREEN</b>	<p>Three-wing hospital screen:</p> <ul style="list-style-type: none"> <li>* Metal construction/ chrome nickel.</li> <li>* Dimensions in unfolded position to 180 cm.</li> <li>* Lightweight construction, to allow active disinfection.</li> <li>* Wheeled base for easy movement and folding of the screen.</li> <li>* The panels shall be made of a material that allows treatment with disinfectants.</li> </ul>	<b>2 Pcs</b>
NO 14 U0066.S		<b>STEP</b>	<p>Medical step:</p> <ul style="list-style-type: none"> <li>* To be of metal tubular construction, powder coated.</li> <li>* Have anti-slip platforms.</li> <li>* The materials of the step to withstand cleaning and disinfection;</li> <li>* Size up to 40/50/42</li> </ul>	<b>1 Pcs</b>
NO 16 U0065.S		<b>EXAMINATION COUCH</b>	<p>Two-section medical couch:</p> <ul style="list-style-type: none"> <li>* Size 180/60/70cm.</li> <li>* Tubular construction of material allowing active disinfection;</li> <li>* To be upholstered in elastic foam and dressed in eco-leather to allow active disinfection</li> <li>* Be able to lift the back section.</li> <li>* Safe Workload: 120 kg.</li> </ul>	<b>3 Pcs</b>
NO 18 U0930		<b>HOSPITAL BED</b>	<p>Stationary hospital bed:</p> <ul style="list-style-type: none"> <li>* Size to 200/90/50 cm.</li> <li>* Complete with passive mattress with size to 200/90/15 cm.</li> </ul>	<b>2 Pcs</b>
NO 19		<b>MATRESSES</b>	<p>Antidekubital mattress made of high quality materials - antibacterial,</p>	<b>2 Pcs</b>

Address:







**HC Medikal A.S.**

06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

[www.ihcmed.com.tr](http://www.ihcmed.com.tr)






# IHC

U0018.S			<p>hypoallergenic, designed for hospital bed;</p> <ul style="list-style-type: none"> <li>* Be watertight;</li> <li>* To be covered with a case, with a flexible polyurethane coating, for maintenance and cleaning with disinfectants;</li> <li>* Overall dimensions up to 200/90/15 cm.</li> </ul>	
NO 20 U0024.S. T8.KL1		<b>BEDSIDE CABINET</b>	<p>Attachment (hospital) locker:</p> <ul style="list-style-type: none"> <li>* The product shall be made of flat material with laminate coating and with stand frequent treatment with disinfectant.</li> <li>* The device includes drawer, door.</li> <li>* Size to 50/50/75 cm.</li> </ul>	<b>2 Pcs</b>
NO 21 U0072		<b>IV STAND</b>	<p>Infusion stand:</p> <ul style="list-style-type: none"> <li>* Metal construction with base on five wheels.</li> <li>* Possibility of mechanical height adjustment;</li> <li>* With four hooks for hanging;</li> <li>* Materials to withstand cleaning and disinfection;</li> </ul>	<b>30 Pcs</b>
NO 22 U0118		<b>PATIENT TRANSPORT CHAIR</b>	<p>Patient transport chair:</p> <ul style="list-style-type: none"> <li>* The product shall be made of stainless steel tubular construction;</li> <li>* The seat and backrest shall be of a polyethylene plastic, for maintenance and cleaning with disinfectants;</li> <li>* With a high garment covering the back and head.</li> <li>* Folding footrests, without protective buffers.</li> <li>* To have comfortable handles for gripping and moving.</li> <li>* To have two bigger back wheels with brake mechanism;</li> <li>* load capacity: 150 kg.</li> </ul>	<b>1 Pcs</b>
NO 23 U0922		<b>CART</b>	<p>Total medicine multipurpose cart- with two levels:</p> <ul style="list-style-type: none"> <li>* To be of metal construction with stainless steel.</li> <li>* To have two stainless steel tops with a bevelled edge;</li> <li>* Wheels to be rubber coated, bearing to allow 360 rotation and chair-pedal.</li> <li>* Materials to withstand frequent treatment with a disinfectant based on sodium precarbonate, quaternary ammonium salts and alkylamines, alcohols and</li> </ul>	<b>2 Pcs</b>





**IHC Medikal A.S.**

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

[www.ihcmed.com.tr](http://www.ihcmed.com.tr)

			didecyldimethylammonium chloride; * Occupying space: to 70/50/90 cm.	
NO 24 U0074		<b>LAUNDRY CART</b>	Laundry cart * To be of sturdy metal construction, pages with metal braiding; * Removable longitudinal half-page have comfortable and easy to move wheels with a load capacity of 250 kg.	<b>2 Pcs</b>
NO 25		<b>TROLLEY</b>	Transport trolley for sanitary and medical materials: Sturdy metal construction made of four sides with metal braid with raster * To have comfortable and easy to move wheels * Load capacity of 150 kg. * Size to 80/60/110 cm.	<b>2 Pcs</b>
NO 28 U0014.S. T4.B2.K0. YWM.W0		<b>HOSPITAL BED WITH 1 CRANK</b>	Two-section hospital bed: * To be with two-section, one mobile and one stationary; * To have steel frame and fundamentals of individual section of steel profile- powder painted construction; * To be able to raise the back section to 70°; * To be installation with soundproof wheels with rubberized coating and stop-system. * Overall dimensions up to 200/100/50 cm.	<b>2 Pcs</b>
NO 29 U0018.S		<b>MATTRESSES</b>	* Antidekubital mattress made of high quality materials - antibacterial, hypoallergenic, designed for hospital bed; * Be watertight; * To be covered with a case, with a flexible polyurethane coating, for maintenance and cleaning with disinfectants; * Overall dimensions up to 200/90/15 cm.	<b>2 Pcs</b>
NO 30 U0024.S. T8.KL1		<b>BEDSIDE CABINET</b>	Attached cabinet: The device is made of material designed for hospital use and can withstand frequent treatment with disinfectants; * The device includes drawer, door. * With wheels, with foot brakes- at two of them;	<b>2 Pcs</b>




# IHC

			* Size to 50/50/75 cm.	
<p><b>NO 31</b></p> <p><b>U0014.S.</b> <b>T4.B2.K0.</b> <b>YWM.W0</b></p>		<p><b>HOSPITAL BED WITH 1 CRANK</b></p>	<p>Two-section hospital bed:</p> <ul style="list-style-type: none"> <li>* To be with two-section, one mobile and one stationary;</li> <li>* To have steel frame and fundamentals of individual section of steel profile- powder painted construction;</li> <li>* To be able to raise the back section to 70°;</li> <li>* To be installation with soundproof wheels with rubberized coating and stop-system.</li> <li>* Overall dimensions up to 200/100/50 cm.</li> </ul>	<p><b>24 Pcs</b></p>
<p><b>NO 32</b></p> <p><b>U0018.S</b></p>		<p><b>MATTRESSES</b></p>	<ul style="list-style-type: none"> <li>* Anti-decubitus mattress, retrieved from high-quality- antibacterial, hypoallergenic, intended for hospital bed;</li> <li>* To be waterproof;</li> <li>* To be wearing case, with elastic polyurethane coating, for maintenance and cleaning with disinfectants;</li> <li>* Overall dimensions up to 200/90/15 cm.</li> </ul>	<p><b>24 Pcs</b></p>
<p><b>NO 33</b></p> <p><b>U0024.S.</b> <b>T8.KL1</b></p>		<p><b>BEDSIDE CABINET</b></p>	<p>Attachment (hospital) locker: The device is made of material designed for hospital use and can withstand frequent treatment with disinfectants;</p> <ul style="list-style-type: none"> <li>* The device includes drawer, door.</li> <li>* With wheels, with foot brakes- at two of them;</li> <li>* Size to 50/50/75 cm.</li> </ul>	<p><b>24 Pcs</b></p>
<p><b>NO 34</b></p> <p><b>U0099.S.</b> <b>T4S.B2.K0</b> <b>.M1.M3.C</b> <b>1.YWM.</b> <b>W8</b></p>		<p><b>HOSPITAL BED</b></p>	<p>Electric hospital bed with adjustable height: Size: min 2100/950/500-820 mm.</p> <ul style="list-style-type: none"> <li>* Steel frame and bases of the individual sections of galvanized steel mesh- powder painting.</li> <li>* ABS trays with two handles, with the possibility of placing a plate.</li> <li>* Lifting the sections by means of electric motor.</li> <li>* Movement of sections by remote control: <ul style="list-style-type: none"> <li>- Lifting at the back 0°-75°</li> <li>- Lifting at the shin 0°-45°</li> <li>- Lifting at the height- 300 mm.</li> </ul> </li> <li>* Inclination of the base: trendelenburg/antitrendelenburg</li> </ul>	<p><b>3 Pcs</b></p>





**IHC Medikal A.S.**

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

[www.ihcmed.com.tr](http://www.ihcmed.com.tr)

			<p>Remote control with buttons for adjusting sections at the back, buttons for adjusting sections at the hip, buttons for simultaneous movement of sections, buttons for adjusting height and buttons for trendelenburg/antitrendelenburg.</p> <ul style="list-style-type: none"> <li>* 4 wheels rotating on 360° with foot brakes- soundproofing.</li> <li>* Stainless steel adjustable/height/infusion stand.</li> <li>* Without side barriers.</li> </ul>	
<p><b>NO 35</b> <b>U0019.S</b></p>		<p><b>MATTRES S</b></p>	<ul style="list-style-type: none"> <li>* Be with a core made of elastic memory foam with one side;</li> <li>* The elasticity of the foam ensures maximum adaptability of the mattress to the anatomical curves of the human body.</li> <li>* Orthopedic properties of the mattress to help keep the spine in the human body in the correct position, ensuring overall relaxation of the body</li> </ul> <p>The mattress shall be of a robust construction not allowing relaxation in lateral directions.</p> <ul style="list-style-type: none"> <li>* The mattress to be single faced.</li> <li>* The mattress cover should provide good ventilation and air circulation of the air to prevent it from retaining moisture and forming an environment for the development of harmful microorganisms.</li> <li>* Occupying space: 82-100/190-200 cm.</li> </ul>	<p><b>3 Pcs</b></p>
<p><b>NO 36</b> <b>U0024.S.</b> <b>T8.KL1</b></p>		<p><b>BEDSIDE CABINET</b></p>	<p>Attachment (hospital) locker: The device is made of material designed for hospital use and can withstand frequent treatment with disinfectants;</p> <ul style="list-style-type: none"> <li>* The device includes drawer, door</li> <li>* With wheels, with foot brakes- at two of them;</li> <li>* Size to 50/50/75 cm.</li> </ul>	<p><b>3 Pcs</b></p>
<p><b>NO 37</b> <b>U0047.S</b></p>		<p><b>MEDICAL CABINET</b></p>	<p>Double medical cabinet:</p> <ul style="list-style-type: none"> <li>* To be used for storage of medical supplies;</li> <li>* Cabinet to be of glazed doors for maximum visibility;</li> <li>* Mechanical parts to ensure very cyclic operation;</li> <li>* To be with adjustable legs,</li> </ul>	<p><b>1 Pcs</b></p>

# IHC







			<p>antibacterial base, neutral to disinfectants with easy to clean surface;</p> <ul style="list-style-type: none"> <li>* The handles shall have a wear-resistant coating;</li> <li>* With secure door locking;</li> <li>* Materials to withstand cleaning with disinfectants.</li> <li>* Dimensions up to 90/45/ 180 cm.</li> </ul>	
<p><b>NO 38</b> <b>U0049.S</b></p>		<p><b>CABINET</b></p>	<p>Double cabinet for medicines</p> <ul style="list-style-type: none"> <li>* The cabinet shall have an upper glazed part divided into three levels and a lower solid part with two levels;</li> <li>* Mechanical parts to ensure very cyclic operation;</li> <li>* To be with adjustable legs, antibacterial base, neutral to disinfectants with easy to clean surface;</li> <li>* The handles shall have a wear-resistant coating;</li> <li>* With secure door locking;</li> <li>* Materials to withstand cleaning with disinfectants</li> <li>* Dimensions up to 90/45/ 180 cm.</li> </ul>	<p><b>1 Pcs</b></p>
<p><b>NO 39</b> <b>U0114.S</b></p>		<p><b>GYNECOLOGICAL CHAIR</b></p>	<p>Gynecological electric chair:</p> <ul style="list-style-type: none"> <li>* Electrostatic powder coating of components:</li> <li>* Electric adjustment of the back section , height trendelenburg;</li> <li>* Mattress platforms, made from polyurethane foam with high density;</li> <li>* Coating of artificial leather, sustainable of fire and disinfectants;</li> <li>* Kit: footrest, headrest, armrests, handles, waste container.</li> <li>* Workload: 180 kg.</li> </ul>	<p><b>1 Pcs</b></p>
<p><b>NO 41</b> <b>U0065.S</b></p>		<p><b>EXAMINATION COUCH</b></p>	<p>Two-section medical couch:</p> <ul style="list-style-type: none"> <li>* Size 180/60/70 cm.</li> <li>* Tubular construction of material allowing active disinfection;</li> <li>* To be upholstered in elastic foam and dressed in eco-leather to allow active disinfection</li> <li>* Be able to lift the back section.</li> <li>* Safe Workload: 120 kg.</li> </ul>	<p><b>2 Pcs</b></p>
<p><b>NO 42</b> <b>U0068.S</b></p>		<p><b>SCREEN</b></p>	<p>Three-wing hospital screen:</p> <ul style="list-style-type: none"> <li>* Metal construction/ chrome nickel .</li> <li>* Dimensions in unfolded position to 180 cm.</li> <li>* Lightweight construction, to allow active disinfection.</li> </ul>	<p><b>2 Pcs</b></p>

**IHC Medikal A.S.**

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

[www.ihcmed.com.tr](http://www.ihcmed.com.tr)





# IHC

			<ul style="list-style-type: none"> <li>* Wheeled base for easy movement and folding of the screen.</li> <li>* The panels shall be made of a material that allows treatment with disinfectants.</li> </ul>	
NO 43 U0066.S		STEP	<p>Medical step:</p> <ul style="list-style-type: none"> <li>* To be of metal tubular construction, powder coated.</li> <li>* Have anti-slip platforms.</li> <li>* The materials of the step to withstand cleaning and disinfections;</li> <li>* Size up to 40/50/42.</li> </ul>	1 Pcs
NO 45 U0065		EXAMINATION COUCH	<p>Two-section medical couch:</p> <p>Size 180/60/70 cm.</p> <ul style="list-style-type: none"> <li>* Tubular construction of material allowing active disinfection;</li> <li>* To be upholstered in elastic foam and dressed in eco-leather to allow active disinfection.</li> <li>* Be able to lift the back section.</li> <li>* Safe Workload: 120 kg.</li> </ul>	2 Pcs
NO 46 U0066.S		STEP	<p>Medical step</p> <ul style="list-style-type: none"> <li>* To be of metal tubular construction, powder coated.</li> <li>* Have anti-slip platforms.</li> <li>* The materials of the step to withstand cleaning and disinfections;</li> <li>* Size up to 40/50/42.</li> </ul>	1 Pcs
NO 48 U0930		HOSPITAL BED	<p>Stationary hospital bed</p> <ul style="list-style-type: none"> <li>* Size to 200/90/50cm.</li> <li>* Complete with passive mattress with size to 200/90/15 cm.</li> </ul>	2 Pcs
NO 49 U0018.S		MATRESSES	<p>Antidekubital mattress made of high quality materials - antibacterial, hypoallergenic, designed for hospital bed;</p> <ul style="list-style-type: none"> <li>* Be watertight;</li> <li>* To be covered with a case, with a flexible polyurethane coating, for maintenance and cleaning with disinfectants;</li> <li>* Overall dimensions up to 200/90/15 cm.</li> </ul>	2 Pcs
NO 50 U0024.S. T8.KL1		BEDSIDE CABINET	<p>Attachment (hospital) locker the device is made of material designed for hospital use and can withstand frequent treatment with disinfectants;</p>	2 Pcs




**IHC Medikal A.S.**

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

[www.ihcmed.com.tr](http://www.ihcmed.com.tr)

			<ul style="list-style-type: none"> <li>* The device includes drawer, door.</li> <li>* With wheels, with foot brakes- at two of them;</li> <li>* Size to 50/50/75 cm.</li> </ul>	
<p><b>NO 51</b></p> <p><b>U0014.S.</b> <b>T4.B2.K0.</b> <b>YWM.WO</b></p>		<p><b>HOSPITAL BED WITH 1 CRANK</b></p>	<p>Two-section hospital bed:</p> <ul style="list-style-type: none"> <li>* To be with two-section, one mobile and one stationary;</li> <li>* To have steel frame and fundamentals of individual section of steel profile- powder painted construction;</li> <li>* To be able to raise the back section to 70°;</li> <li>* To be installation with soundproof wheels with rubberized coating and stop-system.</li> <li>* Overall dimensions up to 210/100/50 cm.</li> </ul>	<p><b>5 Pcs</b></p>
<p><b>NO 52</b></p> <p><b>U0018.S</b></p>		<p><b>MATTRESSES</b></p>	<ul style="list-style-type: none"> <li>* Antidekubital mattress made of high quality materials - antibacterial, hypoallergenic, designed for hospital bed;</li> <li>* Be watertight;</li> <li>* To be covered with a case, with a flexible polyurethane coating, for maintenance and cleaning with disinfectants;</li> <li>* Overall dimensions up to 200/90/15 cm.</li> </ul>	<p><b>5 Pcs</b></p>
<p><b>NO 53</b></p> <p><b>U0024.T8</b> <b>.KL1</b></p>		<p><b>BEDSIDE CABINET</b></p>	<p>Attachment (hospital) locker: the device is made of material designed for hospital use and can withstand frequent treatment with disinfectants;</p> <ul style="list-style-type: none"> <li>* The device includes drawer, door.</li> <li>* With wheels, with foot brakes- at two of them;</li> <li>* Size to 50/50/75 cm.</li> </ul>	<p><b>5 Pcs</b></p>
<p><b>NO 54</b></p> <p><b>U0014.T4</b> <b>.B2.K0.Y</b> <b>WM.WO</b></p>		<p><b>HOSPITAL BED WITH 1 CRANK</b></p>	<p>Two-section hospital bed:</p> <ul style="list-style-type: none"> <li>* To be with two-section, one mobile and one stationary;</li> <li>* To have steel frame and fundamentals of individual section of steel profile- powder painted construction;</li> <li>* To be able to raise the back section to 70°;</li> <li>* To be installation with soundproof wheels with rubberized coating and stop-system.</li> <li>* Overall dimensions up to 210/100/50 cm.</li> </ul>	<p><b>6 Pcs</b></p>





# IHC




<p>NO 55 U0018.S</p>		<p><b>MATTRES S</b></p>	<p>* Antidekubital mattress made of high quality materials - antibacterial, hypoallergenic, designed for hospital bed; * Be watertight; * To be covered with a case, with a flexible polyurethane coating, for maintenance and cleaning with disinfectants; * Overall dimensions up to 200/90/15cm</p>	<p><b>6 Pcs</b></p>
<p>NO 56 U0024.T8 .KL1</p>		<p><b>BEDSIDE CABINET</b></p>	<p>Attachment (hospital) locker The device is made of material designed for hospital use and can withstand frequent treatment with disinfectants; * The device includes drawer, door * With wheels, with foot brakes- at two of them; * Size to 50/50/75 cm.</p>	<p><b>2 Pcs</b></p>
<p>NO 57 U0099.S. T4.B2.K1. M1.M3.C 1.YWM. W8</p>		<p><b>HOSPITAL BED</b></p>	<p>Electric hospital bed with adjustable height: Size: min 2100/950/500-820mm * Steel frame and bases of the individual sections of galvanized steel mesh- powder painting * ABS trays with two handles, with the possibility of placing a plate. * Lifting the sections by means of electric motor * Movement of sections by remote control: - Lifting at the back 0°-75° - Lifting at the shin 0°-45° - Lifting at the height- 300mm. * Inclination of the base: trendelenburg/antitrendelenburg Remote control with buttons for adjusting sections at the back, buttons for adjusting sections at the hip, buttons for simultaneous movement of sections, buttons for adjusting height and buttons for trendelenburg/antitrendelenburg * 4 wheels rotating on 360° with foot brakes- soundproofing * Stainless steel adjustable/height/infusion stand * Side barriers.</p>	<p><b>2 Pcs</b></p>
<p>NO 58 U0019.S</p>		<p><b>MATTRES S</b></p>	<p>* Be with a core made of elastic memory foam with one side; * The elasticity of the foam ensures</p>	<p><b>2 Pcs</b></p>

**IHC Medikal A.S.**






Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40





[www.ihcmed.com.tr](http://www.ihcmed.com.tr)

			<p>maximum adaptability of the mattress to the anatomical curves of the human body</p> <ul style="list-style-type: none"> <li>* Orthopedic properties of the mattress to help keep the spine in the human body in the correct position, ensuring overall relaxation of the body</li> <li>The mattress shall be of a robust construction not allowing relaxation in lateral directions</li> <li>* The mattress to be single faced</li> <li>* The mattress cover should provide good ventilation and air circulation of the air to prevent it from retaining moisture and forming an environment for the development of harmful microorganisms.</li> <li>* Occupying space: 82-100/190-200 cm.</li> </ul>	
<p><b>NO 59</b></p> <p><b>U0024.T8</b></p> <p><b>.KL1</b></p>		<p><b>BEDSIDE CABINET</b></p>	<p>Attachment (hospital) locker</p> <p>The device is made of material designed for hospital use and can withstand frequent treatment with disinfectants;</p> <ul style="list-style-type: none"> <li>* The device includes drawer, door</li> <li>* With wheels, with foot brakes- at two of them;</li> <li>* Size to 50/50/75cm.</li> </ul>	<p><b>2 Pcs</b></p>
<p><b>NO 60</b></p> <p><b>U0014.T4</b></p> <p><b>.B2.K0.Y</b></p> <p><b>WM.WO</b></p>		<p><b>HOSPITAL BED WITH 1 CRANK</b></p>	<p>Two-section hospital bed:</p> <ul style="list-style-type: none"> <li>* To be with two-section, one mobile and one stationary;</li> <li>* To have steel frame and fundamentals of individual section of steel profile- powder painted construction;</li> <li>* To be able to raise the back section to 70°;</li> <li>* To be installation with soundproof wheels with rubberized coating and stop-system.</li> <li>* Overall dimensions up to 210/100/50 cm.</li> </ul>	<p><b>2 Pcs</b></p>
<p><b>NO 61</b></p> <p><b>U0018.S</b></p>		<p><b>MATTRES S</b></p>	<p>Antidekubital mattress made of high quality materials - antibacterial, hypoallergenic, designed for hospital bed;</p> <ul style="list-style-type: none"> <li>* Be watertight;</li> <li>* To be covered with a case, with a flexible polyurethane coating, for maintenance and cleaning with</li> </ul>	<p><b>2 Pcs</b></p>


			<p>disinfectants;          * Overall dimensions up to 200/90/15cm.</p>	
<p><b>NO 62</b>  <b>U0024.T8</b>  <b>.KL1</b></p>		<p><b>BEDSIDE CABINET</b></p>	<p>Attachment (hospital) locker:          The device is made of material designed for hospital use and can withstand frequent treatment with disinfectants;          * The device includes drawer, door          * With wheels, with foot brakes- at two of them;          * Size to 50/50/75cm.</p>	<p><b>2 Pcs</b></p>
<p><b>NO 63</b>  <b>U0047</b></p>		<p><b>MEDICAL CABINET</b></p>	<p>Double medical cabinet          * To be used for storage of medical supplies;          * Cabinet to be of glazed doors for maximum visibility;          * Mechanical parts to ensure very cyclic operation;          * To be with adjustable legs, antibacterial base, neutral to disinfectants with easy to clean surface;          * The handles shall have a wear-resistant coating;          * With secure door locking;          * Materials to withstand cleaning with disinfectants          * Dimensions up to 90/45/ 180cm</p>	<p><b>1 Pcs</b></p>
<p><b>NO 64</b>  <b>U0049</b></p>		<p><b>CABINET</b></p>	<p>Double medicine cabinet          * The cabinet shall have an upper glazed part divided into three levels and a lower solid part with two levels;          * Mechanical parts to ensure very cyclic operation;          * To be with adjustable legs, antibacterial base, neutral to disinfectants with easy to clean surface;          * The handles shall have a wear-resistant coating;          * With secure door locking;          * Materials to withstand cleaning with disinfectants          * Dimensions up to 90/45/ 180cm</p>	<p><b>1 Pcs</b></p>

<p>NO 65 U0114.S</p>		<p><b>GYNECOLOGICAL CHAIR</b></p>	<ul style="list-style-type: none"> <li>* Electrostatic powder coating of components;</li> <li>* Electric adjustment of the back section, height, trendelenburg;</li> <li>* Mattress platforms, made from polyurethane foam with high density;</li> <li>* Coating of artificial leather, sustainable of fire and disinfectants;</li> <li>* Kit: footrest, headrest, armrests, handles, waste container</li> <li>* Workload: 180 kg.</li> </ul>	<p><b>1 Pcs</b></p>
<p>NO 67 U0065</p>		<p><b>EXAMINATION COUCH</b></p>	<p>Two-section medical couch: Size 180/60/70cm.</p> <ul style="list-style-type: none"> <li>* Tubular construction of material allowing active disinfection;</li> <li>* To be upholstered in elastic foam and dressed in eco-leather to allow active disinfection</li> <li>* Be able to lift the back section.</li> <li>* Safe Workload: 120 kg.</li> </ul>	<p><b>1 Pcs</b></p>
<p>NO 68 U0068</p>		<p><b>SCREEN</b></p>	<p>Three-wing hospital screen:</p> <ul style="list-style-type: none"> <li>* Metal construction/ chrome nickel.</li> <li>* Dimensions in unfolded position to 180cm.</li> <li>* Lightweight construction, to allow active disinfection.</li> <li>* Wheeled base for easy movement and folding of the screen.</li> <li>* The panels shall be made of a material that allows treatment with disinfectants.</li> </ul>	<p><b>1 Pcs</b></p>
<p>NO 69 U0066</p>		<p><b>STEP</b></p>	<p>Medical step:</p> <ul style="list-style-type: none"> <li>* To be of metal tubular construction, powder coated.</li> <li>* Have anti-slip platforms.</li> <li>* The materials of the step to withstand cleaning and disinfection;</li> <li>* Size up to 40/50/42</li> </ul>	<p><b>1 Pcs</b></p>
<p>NO 71 U0049</p>		<p><b>MEDICAL CABINET</b></p>	<p>Double medicine cabinet:</p> <ul style="list-style-type: none"> <li>* The cabinet shall have an upper glazed part divided into three levels and a lower solid part with two levels;</li> <li>* Mechanical parts to ensure very cyclic operation;</li> <li>* To be with adjustable legs, antibacterial base, neutral to disinfectants with easy to clean surface;</li> <li>* The handles shall have a wear-resistant coating;</li> </ul>	<p><b>1 Pcs</b></p>

			<ul style="list-style-type: none"> <li>* With secure door locking;</li> <li>* Materials to withstand cleaning with disinfectants</li> <li>* Dimensions up to 90/45/ 180cm</li> </ul>	
<p><b>NO 72</b></p> <p><b>U0114.S</b></p>		<p><b>GYNECOLOGICAL CHAIR</b></p>	<p>Gynecological mechanic chair:</p> <ul style="list-style-type: none"> <li>* Electrostatic powder coating of components:</li> <li>* Electric adjustment of the back section, height, trendelenburg;</li> <li>* Mattress platforms, made from polyurethane foam with high density;</li> <li>* Coating of artificial leather, sustainable of fire and disinfectants;</li> <li>* Kit: footrest, headrest, armrests, handles, waste container</li> <li>* Workload: 180 kg.</li> </ul>	<p><b>1 Pcs</b></p>
<p><b>NO 74</b></p> <p><b>U0065.S</b></p>		<p><b>EXAMINATION COUCH</b></p>	<p>Two-section medical couch: Size 180/60/70cm.</p> <ul style="list-style-type: none"> <li>* Tubular construction of material allowing active disinfection;</li> <li>* To be upholstered in elastic foam and dressed in eco-leather to allow active disinfection</li> <li>* Be able to lift the back section.</li> <li>* Safe Workload: 120 kg.</li> </ul>	<p><b>1 Pcs</b></p>
<p><b>NO 75</b></p> <p><b>U0066.S</b></p>		<p><b>STEP</b></p>	<p>Medical step:</p> <ul style="list-style-type: none"> <li>* To be of metal tubular construction, powder coated.</li> <li>* Have anti-slip platforms.</li> <li>* The materials of the step to withstand cleaning and disinfection;</li> <li>* Size up to 40/50/42.</li> </ul>	<p><b>1 Pcs</b></p>
<p><b>NO 77</b></p> <p><b>U0930</b></p>		<p><b>HOSPITAL BED</b></p>	<p>Stationary hospital bed:</p> <ul style="list-style-type: none"> <li>* Size to 200/90/50cm.</li> <li>* Complete with passive mattress with size to 200/90/15cm.</li> </ul>	<p><b>1 Pcs</b></p>
<p><b>NO 78</b></p> <p><b>U0018.S</b></p>		<p><b>MATTRESSES</b></p>	<p>Antidekubital mattress made of high quality materials - antibacterial, hypoallergenic, designed for hospital bed;</p> <ul style="list-style-type: none"> <li>* Be watertight;</li> <li>* To be covered with a case, with a flexible polyurethane coating, for maintenance and cleaning with disinfectants;</li> <li>* Overall dimensions up to 200/90/15cm</li> </ul>	<p><b>1 Pcs</b></p>

<p><b>NO 79</b> <b>U0024.T8</b> <b>.KL1</b></p>		<p><b>BEDSIDE CABINET</b></p>	<p>Attachment (hospital) locker: The device is made of material designed for hospital use and can withstand frequent treatment with disinfectants; * The device includes drawer, door * With wheels, with foot brakes- at two of them; * Size to 50/50/75cm</p>	<p><b>1 Pcs</b></p>
<p><b>NO 80</b> <b>U0072</b></p>		<p><b>IV STAND</b></p>	<p>Infusion stand: * Metal construction with base on five wheels. * Possibility of mechanical height adjustment; * with four hooks for hanging; * materials to withstand cleaning and disinfections;</p>	<p><b>40 Pcs</b></p>
<p><b>NO 81</b> <b>U0118.S</b></p>		<p><b>PATIENT TRANSPORT CHAIR</b></p>	<p>Patient transport chair: * The product shall be made of stainless steel tubular construction; * The seat and backrest shall be of a polyethylene plastic, for maintenance and cleaning with disinfectants; * With a high garment covering the back and head. * Folding footrests, without protective buffers. * To have comfortable handles for gripping and moving. * To have two bigger back wheels with brake mechanism; * load capacity: 150 kg.</p>	<p><b>1 Pcs</b></p>
<p><b>NO 82</b> <b>U0922</b></p>		<p><b>CART</b></p>	<p>Total medicine multipurpose cart- with two levels: * To be of metal construction with chrome plating. * To have two stainless steel tops with a bevelled edge; * Wheels to be rubber coated, bearing to allow 360 rotation and chair-pedal. * Materials to withstand frequent treatment with a disinfectant based on sodium precarbonate, quaternary ammonium salts and alkylamines, alcohols and didecylmethylammonium chloride; * Occupying space: to 70/50/90cm.</p>	<p><b>2 Pcs</b></p>

# IHC

<p><b>NO 83</b></p> <p><b>U0074</b></p>		<p><b>LAUNDRY CART</b></p>	<p>Laundry cart</p> <ul style="list-style-type: none"><li>* To be of sturdy metal construction, pages with metal braiding;</li><li>* Removable longitudinal half-page have comfortable and easy to move wheels with a load capacity of 250 kg</li></ul>	<p><b>2 Pcs</b></p>
---	---	--------------------------------	---	---------------------

**IHC Medikal A.S.**

Address: Bahcelievler District Yavuz Road No: 31 06850 Hasanoglan/Elmadag/Ankara/Turkey - Phone: +90 312 866 40 40

[www.ihcmed.com.tr](http://www.ihcmed.com.tr)