


To: UMBAL D-r Georgi Stranski


Pleven, Republic of Bulgaria

Re: announcement No. 00178-2023-0024

for public procurement procedure

## CATALOG

<p>U0973.T4. BS.K4.YS. W8S</p>		<p><b>HOSPITAL BED WITH 1 CRANK</b></p>	<p><i>Mechanical hospital bed with one function.</i></p> <ul style="list-style-type: none"><li>• <i>Mechanical with one function.</i></li><li>• <i>Overall Dimensions: 2050/860/450mm +/-1-2 mm/</i></li><li>• <i>Dimensions of the inner frame: 1990/800 mm +/-1-2 mm/.</i></li><li>• <i>Maximum safe load: 120kg.</i></li><li>• <i>Description:</i></li><li>• <i>Two-section internal frame made of metal CB profiles - powder coating.</i></li><li>• <i>External supporting frame of CB profiles with dimensions of min. 60x30mm and internal frame of transverse CB profiles with dimensions of min. 30x20mm.</i></li><li>• <i>Legs from CB profile with dimensions min 50x50mm.</i></li><li>• <i>Metal frame made of a tube with a diameter of min. f33.7mm at the head and legs welded to the main structure - non-removable for greater strength, with chipboard boards - possibility to choose the color of the chipboard boards.</i></li><li>• <i>Lifting the section by means of 1pc. manual screw mechanism with chrome handle that can be stored under the footboard.</i></li><li>• <i>Movement of the sections: lifting at the back 0°-75°.</i></li><li>• <i>4 pcs. wheels f125mm +/-1-5 mm/., with rubber coating, rotating 360°, each with a foot brake - isolated with a plastic fork, preventing contamination.</i></li><li>• <i>Aluminum side partitions with a lifting and lowering mechanism with dimensions in the unfolded position min. 380mm in height and min. 1470mm. in length.</i></li><li>• <i>Adjustable/in height from 890 to 1590mm//+/-1-5 mm/. infusion stand with four hooks, made of stainless steel, with the possibility of being placed on all four sides of the bed.</i></li><li>• <i>Possibility to place the lifting stand on both sides at the head of the hospital bed.</i></li></ul>	<p><b>25 Pc</b></p>
--	---	---	---	-------------------------

<p><b>U0018.S</b></p>		<p><b>ANTI DECUBITUS MATTRESS</b></p>	<ul style="list-style-type: none"> <li>• <i>Dimensions: 1950x800x120mm /corresponding to the size of the bed/.</i></li> <li>• <i>Interior:</i></li> <li>• <i>foam - anti-decubitus, cut into cubes with a size of min. 70/70mm.</i></li> <li>• <i>foam density min. T25.</i></li> <li>• <i>Cover fabric:</i></li> <li>• <i>Double elastic polyurethane coating.</i></li> <li>• <i>The fabric must meet the highest requirements for optimal weight and resistance to the smallest particles.</i></li> <li>• <i>Waterproof and "breathable" fabric (vapor permeable).</i></li> <li>• <i>Antibacterial, hypoallergenic.</i></li> <li>• <i>Washing at a temperature of 95°C.</i></li> <li>• <i>To withstand modern disinfectants.</i></li> <li>• <i>With zipper (on the short and long sides).</i></li> <li>• <i>Ecological.</i></li> <li>• <i>CRIB 7 fire resistant.</i></li> </ul>	<p><b>25 Pc</b></p>
<p><b>W2S</b></p>		<p><b>LIFTING STAND WITH ADJUSTABLE HANDLE</b></p>	<ul style="list-style-type: none"> <li>• <i>Metal construction - powder coated.</i></li> <li>• <i>Plastic handle with the possibility of adjusting the length by means of a button - with a minimum stroke of 300 mm.</i></li> <li>• <i>The possibility of placing the stand on both sides of the hospital bed and the possibility of 360 degree rotation.</i></li> </ul>	<p><b>3 Pc</b></p>

**U0003.S.T1  
6S.K5.E1.M  
1.M3.C1.Y  
WMS.W1**




**HOSPITAL BED  
WITH 4  
MOTORS**

*Resuscitation medical hospital bed with ABS partitions - four-motor.*

- *Electric with five functions.*
- *Dimensions: 2120/1000/500-820 mm.*
- *Internal frame dimensions: 2000/840 mm.*
- *Maximum safe load: 200 kg*
- *Description:*
- *Four-section internal frame - one stationary section and three movable*
- *Steel supporting frame and foundations of the individual sections made of galvanized steel mesh and metal CB profiles - powder coating according to RAL 9002.*
- *ABS boards with two handles each - removable /removable/, with the possibility of placing a plate.*
- *Lifting of the sections by means of electric motors.*
- *Movement of the sections by remote control.*
- *Back lift 0°-75°.*
- *Lower leg/calf raise 0°-40°.*
- *Lifting height – 320 mm.*
- *Base slope: Trendelenburg/anti-trendelenburg min. ±12°.*
- *Electric motors – 4 pcs.*
- *Push/pull force: min. 6000/3000 N.*
- *Speed: min. 3.5 mm/s.*
- *Electronic block – 1 pc. with a battery for one-time storage of the motor at the back in the event of a power failure.*
- *Remote control with a key to lock the functions if necessary:*
- *2 pcs. back section adjustment buttons.*
- *2 pcs. lower leg section adjustment buttons.*
- *2 pcs. buttons for simultaneous movement of the back and leg sections /autocracy/.*
- *2 pcs. height adjustment buttons.*
- *2 pcs. trendelenburg/anti-trendelenburg buttons.*
- *4 pcs. wheels f 125 mm., with rubber coating, rotating at 3600 each with a foot brake - isolated with a plastic fork, preventing pollution.*
- *ABS side partitions with a mechanism for lowering and raising by means of a shock absorber - 4 pcs.*
- *Adjustable /in height from 890 to 1590 mm/ stainless steel infusion stand with possibility to be placed on both sides of the bed - four hooks.*

**43  
Pc**

<b>U0018.S</b>		<b>ANTI DECUBITUS MATTRESS</b>	<ul style="list-style-type: none"><li>• <i>Dimensions: 2000x880x120 mm /corresponding to the size of the bed/.</i></li><li>• <i>Interior:</i></li><li>• <i>foam - anti-decubitus, cut into cubes of size 70x70.</i></li><li>• <i>foam density min. T25.</i></li><li>• <i>Cover fabric:</i></li><li>• <i>Double elastic polyurethane coating.</i></li><li>• <i>The fabric must meet the highest requirements for optimal weight and resistance to the smallest particles.</i></li><li>• <i>Waterproof and "breathable" fabric (vapor permeable).</i></li><li>• <i>Antibacterial, hypoallergenic.</i></li><li>• <i>Washing at a temperature of 95°C.</i></li><li>• <i>To withstand modern disinfectants.</i></li><li>• <i>With a zipper (on the short and long side).</i></li><li>• <i>Ecological.</i></li><li>• <i>CRIB 7 fire resistant.</i></li></ul>	<b>43 Pc</b>
----------------	---	--	--	------------------