




PRODUCT CODE	PRODUCT PHOTO	PRODUCT NAME	DESCRIPTION
<p>U0014.S.T 36.B2.KS.Y WM.W8S</p> <p>U0018.S</p>	 	<p>HOSPITAL BED WITH 1 CRANK</p> <p>ANTI DECUBITUS MATTRESS</p>	<p>I Mechanical hospital bed with sideboards, infusion stand and mattress.</p> <p><u>1 Hospital bed:</u></p> <p>1.1 Overall Dimensions: 2120/970/450 mm. 1.2 Dimensions of internal frame: 2000/900 mm. 1.3 Maximum safe load: 150 kg. 1.4 Description:</p> <ul style="list-style-type: none"> * Carrying frame made of steel profiles. * Two-section internal frame made of metal profiles. * Electrostatic powder painting. - Possibility to choose a color according to the RAL system. * Manual mechanism for lifting the section at the back. - Screw. - Chrome handle with the possibility of being stored under the bed frame. * Movement of the back section: 70°. * Plastic boards: - Removable for easy access to the patient. - With two handles for easy positioning/moving of the bed. - Possibility to place a plate at the footboard. - Possibility to choose a color / blue or beige decoration /. * Side partitions: - Aluminum. - Lifting and lowering mechanism. * Infusion stand: - Adjustable in height from 890 to 1590 mm. - Material: stainless steel. - Four hooks. - Possibility of placement on both sides of the bed. *Wheels: - 4 pcs. /removable, mounted with a bolt to the metal structure, allowing easy replacement in case of damage/. - Size: 125 mm. - 360° rotating. - Foot brake on each wheel. - Rubber coating. - Metal fork with plastic coating on the wheel, preventing pollution. <p><u>2 Hospital mattress:</u></p> <p>2.1 Dimensions: 2000x880x120 mm. 2.2 Interior:</p> <ul style="list-style-type: none"> - Material: polyurethane foam. - Density 25 kg/m3. - Hardness 4.2 kPa. - Antidecubitus - cut into cubes of size 70/70 mm. <p>2.3 Cover Fabric:</p>

			<ul style="list-style-type: none"> - Double elastic polyurethane coating on a polyester base. - The fabric meets the highest requirements for optimal weight and resistance to the smallest particles. - Waterproof and breathable fabric (vapor permeable). - Water vapor permeability: min. < 1950 g/m² - 24 h at 38° C. - Antibacterial, hypoallergenic. - Washing at a temperature of 95°C. - Resists modern disinfectants containing. - Denatured alcohol: 90% alcohol + 2% ether. - With a zipper (on the short and long side). - CRIB 7 fire rated - BS 7175 standard. - Ability to choose the color of the case.
<p>U0112.S.T 8.KL2</p>		<p>BEDSIDE CABINET</p>	<p>1 Dimensions: 465/465/795 mm. 2 Beige color. 3 Material: ABS plastic, the corners aluminium. 4 Departments: - A cabinet with a compartment divided by a shelf, the possibility of placing the shelf in at least 2 positions. - A door with a compartment for belongings on the inside. - A drawer. - Movable tray under the top of the cabinet. 5 Towel hangers on the left and right side of the locker - 4 pcs. 6 Wheels: - 4 pcs, 50 mm twin, plastic wheels. - 360° rotating. - Removable/removable/- allowing easy replacement in case of damage. - With foot brakes on two of them.</p>

IHC MEDİKAL A.Ş.
Hasanoğlu Bahçelievler Mh. Yavuz Cad.
No:31 Elmadağ/ANKARA
Gsm : 0544 414 72 42
Elmadağ V.D. 465 127 1713