

U0069

SCREEN

GENERAL SPECIFICATIONS

- Four winged screen is adjusted by castors.
- The curtain is removable, attachable and washable.
- Epoxy painted steel construction is resistant to corrosion.
- Screen is able to be used in narrow area, owing to foldable wings system if necessary.

TECHNICAL SPECIFICATIONS

- | | |
|----------|-----------------------|
| ■ Length | : 280 cm |
| ■ Height | : 166 cm |
| ■ Weight | : approx 13,8 kg |
| ■ Volume | : 0,15 m ³ |



U0070

WHEELCHAIR

GENERAL SPECIFICATIONS

- With foldable property, wheelchair is able to be placed in narrow area.
- Moving foot place.
- Epoxy painted steel construction is resistant to corrosion.
- With an effective brake system, able to stopped easily.
- Back and front wheels have filling properties.

TECHNICAL SPECIFICATIONS

- | | |
|------------------|-----------------------|
| ■ Width | : 62,5 cm |
| ■ Depth | : 107 cm |
| ■ Sitting Height | : 50 cm |
| ■ Height | : 94 cm |
| ■ Weight | : approx 22 kg |
| ■ Volume | : 0,21 m ³ |

