

U0075

LINEN TROLLEY CLEAN & DIRTY

GENERAL SPECIFICATIONS

- Made of heat and water resistant compact wood housing made.
- Three shelves for clean laundry.
- Can be collapsible, washable laundry bag.
- Made of aluminum that can be easily switched on and off dirty laundry cover system.
- There are bumper on corners to protect.
- Moves easily by 100 mm diameter four castors. Two of them with brake to fix the position.

TECHNICAL SPECIFICATIONS

- | | |
|----------|-----------------------|
| ■ Length | : 95 cm |
| ■ Depth | : 66,5 cm |
| ■ Height | : 107 cm |
| ■ Weight | : approx 55 kg |
| ■ Volume | : 0,73 m ³ |

