

GENERAL SPECIFICATIONS

- Resuscitation trolley for emergency.
- mild steel construction.
- ABS plastic top.
- 4 drawers.
- Overall Dimensions: 53 X 80 X 95 cm.
- 4 drawers with central lock.
- Ergonomic Push handle at 1 side.
- 125 mm antistatic castors with locking.
- Cardiac board housing brackets at rear of trolley.
- Universal rail system fitted to width of trolley.
- Pull-out Shelf.
- Overhead frame for some of the accessories.
- Standard Accessories :
  - 1 x Aneroid sphygmomanometer.
  - 1 x Side Accessory mount rail.
  - 1 x Cylinder holder bracket, 48 cft.
  - 1 x Double hook stainless steel I.V. Pole.
  - 1 x Cardiac board.
  - 1 x Swivel Shelf (Monitor/Defibrillator).
  - 1 x Waste Basket/Sharp Container.
  - 1 x Drawer Dividers.



Images are for illustrative purposes only.  
The actual product may differ based on the configuration.