U0012 HOSPITAL BED WITH 1 CRANK

GENERAL SPECIFICATIONS

- Easily cleanable metal bed platform designed for the ventilation of the bed.
- Head movement with rotating levers.
- Ergonomic, durable and easily cleanable metal rails system with maximum safety.
- Melamine coated chipboard head and foot panels.
- Metal components dyed with anti-corrosion electrostatic powder paint after being treated with phosphate.
- Optional melamine coated chipboard color choice.(RSX)

TECHNICAL SPECIFICATIONS

Width	:	102 cm
Length	:	208 cm
Height	:	50 cm
Backrest Adjustment	:	70°
Weight	:	~ 72 kg
Volume	:	0,64 m³

OPTIONS

GROUP 1 / CASTOR OPTIONS



GROUP 2 / SIDE RAIL OPTIONS



GROUP 3 / OTHER ACCESSORY OPTIONS



DESCRIPTION OF OPTION CODES

- T0 WITHOUT CASTORS
- T4 METAL CASTOR TWO WITH BRAKE (125 MM)
- T16 CASTOR TWO WITH BRAKE (125 MM)
- K0 WITHOUT SIDE RAIL
- K1 WITH METAL SIDE RAIL
- K3 WITH METAL SIDE RAIL
- W0 WITHOUT OTHER ACCESSORIES
- W2 LIFTING POLE : BEARING CAPACITY 75 Kg
- W8 IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM

Must select only one of them.

Do not select the others if you select without accessory.
If you select the others you can select all or 1-2-3 ... of them.



* RS : Melamine color options page 83.