## U0018 MATTRESS

### GENERAL SPECIFICATIONS

- High density, specially cut, sliced sponge that can be easily formed in accordance with the bed functions.
- Anti-bacterial bed cover provides air circulation and protects the bed with its liquid tight feature.
- Bed cover, which can be removable and assembled by a zip fastener system if necessary, is of washable feature.

### TECHNICAL SPECIFICATIONS

| | | |
|---|---|---|
| Width | : | 85 cm |
| Depth | : | 194 cm |
| Height | : | 12 cm |
| Weight | : | ~ 8 kg |
| Volume | : | 0,21 $m^3$ |



## U0019 MATTRESS VISCO

### GENERAL SPECIFICATIONS

- Special foam sponge is able to be take a shape according to the hospital bed positions.
- Anti-bacterial bed cover provides air circulation and protects the bed with its liquid tight feature.
- Bed cover, which can be removable and assembled by a zip fastener system if necessary, is of washable feature.

### TECHNICAL SPECIFICATIONS

| | | |
|---|---|---|
| Width | : | 85 cm |
| Depth | : | 194 cm |
| Height | : | 12 cm |
| Weight | : | ~ 8 kg |
| Volume | : | 0,21 $m^3$ |

