U0056 BLOOD RECEPTION CHAIR

GENERAL SPECIFICATIONS

- Antibacterial and cleanable artificial leather covered high density foam cover.
- Height and back part chair adjustable by gas spring.
- Height and angle of arm holder adjustable by mechanical system.
- Epoxy painted steel construction is resistant to corrosion.
- Moves easily by 50 mm diameter castors.
- Vinlex colors are optional.(RVX)

TECHNICAL SPECIFICATIONS

Width	:	43 cm
Depth	:	48 cm
Sitting Height	:	60 cm
Weight	:	approx 13 kg
Volume	:	0,23 m^3



* RV : Vinlex color options page 83.

U0057 BLOOD RECEPTION CHAIR

GENERAL SPECIFICATIONS

- Antibacterial and cleanable artificial leather covered high density foam cover.
- Backwards at the same time moving the foat platform.
- Height and angle of arm holder adjustable by mechanical system.
- Height adjustable IV stand.
- Epoxy painted steel construction is resistant to corrosion.
- Optional artifical leather choice. (RDX)

TECHNICAL SPECIFICATIONS

Width	:	70,5 cm
Depth	:	96 - 170 cm
Sitting Height	:	44 cm
Backrest Adjustment	:	25° - 65°
Weight	:	approx 33 kg
Volume	:	0,90 m^3



* RD : Artifical leather color options page 84.

ISO 9001:2015 ISO 14001:2015 CE