U0108.S MEDICINE TROLLEY

GENERAL SPECIFICATIONS

- Medicine Trolley
- For the storage of medication by ward nursing staff.
- Different compartments for medicines with labelling.
- Body structure made from Mild Steel powder coated.
- ABS plastic top.
- Central locking for all drawers.
- Lockable Castors.
- Rubber Bumpers all around.
- Size: approx. 26" L x 20" W x 40" H.
- Side folding shelf.
- Mobile on 125mm casters, diagonal two with brakes.



Images are for illustrative purposes only.
The actual product may differ based on the configuration.